



## OPTIMIZATION OF THE SEPARATION OF TOCOPHEROLS USING SUB-ZERO TEMPERATURE

### INTRODUCTION

Vitamin E is a term used to describe compounds that contain the 6-hydroxychroman ring and possess the biological activity of  $\alpha$ -tocopherol. There are 8 naturally occurring homologs comprised of 4 tocopherols ( $\alpha$ -,  $\beta$ -,  $\delta$ -,  $\gamma$ -) and 4 tocotrienols ( $\alpha$ -,  $\beta$ -,  $\delta$ -,  $\gamma$ -). In addition to their vitamin E activity, these compounds have been implicated in the reduced risk of many diseases including cancer, cardiovascular disease, Alzheimer's, and cataracts. Vitamin E also provides the greatest total antioxidant activity in the lipoproteins in blood.

While  $\alpha$ -tocopherol has the greatest amount of vitamin E activity, the other homologs may be more effective at specific activities. For example,  $\gamma$ -tocopherol is thought to be more effective at reducing the risk of prostate cancer.

### EXPERIMENTAL

HPLC conditions are summarized in Table 1. One separation was done at ambient temperature and one was performed at  $-20^{\circ}\text{C}$ .

### RESULTS

No separations of all 8 vitamin E homologs have been previously reported using isocratic reversed-phase HPLC. Figure 1 (see next page) illustrates an injection of all 8 homologs on a C18 column using an ambient isocratic separation. Note the coelution of both  $\beta$ - and  $\gamma$ -tocotrienol and of  $\beta$ - and  $\gamma$ -tocopherol. In Figure 2, the separation is performed at  $-20^{\circ}\text{C}$  using the Polaratherm. The low temperature influences the column selectivity resulting in partial separation of the two tocotrienols and near-baseline separation ( $\alpha=1.065$ ) of the two tocopherols.

Table 1:  
Conditions for Analysis of Tocopherols

Column:	Jones Chromatography Genesis C18, 150 x 4.6 mm, 4 $\mu\text{m}$
Mobile Phase:	92:8 acetonitrile:water (ambient temperature) 90:10 acetonitrile:THF ( $-20^{\circ}\text{C}$ )
Flow Rate:	1.0 mL/min
Detection:	Absorbance at 296 nm

### CONCLUSIONS

An isocratic separation of tocopherols at  $-20^{\circ}\text{C}$  was achieved. The separation showed improved resolution of several isomers when compared to the analysis performed at ambient temperature.

### ACKNOWLEDGEMENT

Many thanks to Dr. Neal Craft and John Estes of Craft Technologies, Inc. in Wilson, NC, who developed this separation.



Figure 1: Separation of tocopherols at ambient temperature. Several of the isomers coelute.





Figure 1: Separation of tocopherols at -20 $^{\circ}$ C.  $\beta$ -tocopherol and  $\gamma$ -tocopherol are now well-resolved, and  $\beta$ -tocotrienol and  $\gamma$ -tocotrienol are partially resolved.