

# **Extended Temperature Range Impact on HPLC Column Stability and Performance**

**Stephanie J. Marin**

**Dale Felix, Jody Clark**

**Selerity Technologies, Inc.**

**Salt Lake City, UT**

**[www.selerity.com](http://www.selerity.com)**



# The Selerity Series 9000 Total Temperature Controller

- Forced air oven and chiller
- Isothermal and thermal gradient operation
  - Sub-zero to 200°C
  - Flow rates up to 10.0 mL/min
  - Thermal gradient up to 30°C/min
- Mobile phase pre-heating
- Peltier effluent cooling
- Vapor sensor
- Compatible with any HPLC system



# The High Temperature HPLC Advantage

- **Temperature Programming**
- **Speed**
- **Efficiency**



# Better Chromatography with Temperature Gradient Programming

- **Change retention through temperature gradient programming**
  - Replace solvent gradients with temperature gradients
  - Water less polar and more like methanol so less organic modifier needed



# Faster and More Efficient Separations

- **Speed**
  - Flatter van Deemter curves allow operation at flow rates many times optimal velocity
- **Higher efficiency - better resolution**
  - Increased diffusion rates provide lower plate heights at higher flow rates
  - Lower viscosity and back pressure permits higher flow rates with smaller particle size packings



# The Real Power of Temperature in HPLC is Temperature Programming – but what about columns?

- Evaluated several different types of columns under programmed HTLC conditions
- Compared thermal gradient and solvent gradient results




# Column Characteristics

- Selerity **Blaze** C<sub>8</sub>
  - Polydentate silica, 3 $\mu$ m, 100 x 4.6 mm
  - Stable to 100°C, pH 2-8
- Hamilton PRP<sup>®</sup>-1
  - PSDVB, 5 $\mu$ m, 100 x 4.1 mm
  - Stable to 150°C, pH 0-14
- Thermo Hypersil-Keystone Hypercarb<sup>®</sup>
  - Graphitic carbon, 7 $\mu$ m, 100 x 4.6 mm
  - Stable to 200°C, pH 0-14
- Zirchrom PBD, Diamondbond and CARB



# Traditional Silica Columns Can't Take the Heat


Water attacks siloxane bond or behind point of  
phase attachment






# The Selerity **Blaze** Silica Column Can Take the Heat

- Selerity polydentate phase protects the silanol groups
- Polymer attaches to backbone at several points



# Blank Thermal Gradients

## 50:50 Acetonitrile:Water at 254 nm



# Column Evaluation Conditions

- Seven component mix of acidic, basic and neutral compounds:

**acetophenone**

**amitriptylene**

**aniline**

**Ibuprofen**

**salicylic acid**

**styrene glycol**

**2-phenyl-2-propanol**



# Column Evaluation Conditions

- **Three mobile phases:**
  - acetonitrile:water
  - acetonitrile:water with 0.1% TFA pH  $\approx$  2
  - acetonitrile:ammonium hydroxide pH 10
- **Analyzed using a solvent gradient and a thermal gradient to give similar retention times.**
- **50:50 acetonitrile:aqueous**
  - 50% to 100% ACN over 10 minutes, hold five minutes
  - 50°C to 100°C, 150°C or 200°C at 15°/min, hold five minutes




# Column Evaluation Conditions


- **retention times collected individually**
  - Eight solvent gradients and eight temperature gradients per column
- **column was evaluated initially and between each mobile phase change with two component test mix**
  - uracil and phenol
  - 40°C
  - 50:50 ACN:Water
- **retention time, peak area and efficiency were recorded**




# RT vs °C and RT vs %ACN for PRP<sup>®</sup>-1 ACN:Water



# RT vs °C and RT vs %ACN for PRP®-1 ACN:Water with 0.1% TFA




# RT vs °C and RT vs %ACN for PRP®-1 ACN:Water pH 10







# RT vs °C and RT vs %ACN for Hypercarb<sup>®</sup> Column ACN:Water with 0.1%TFA




# RT vs °C and RT vs %ACN for Hypercarb<sup>®</sup> Column ACN:water pH 10



# RT vs °C and RT vs %ACN for Blaze C<sub>8</sub> ACN:Water



# RT vs °C and RT vs %ACN for Blaze C<sub>8</sub> ACN:TFA




# Peak Capacity Comparison for Thermal and Solvent Gradients

| Mobile Phase                                       | Blaze™ C <sub>8</sub> |                  | PRP-1®           |                  | HyperCarb®       |                  |
|--|-----------------------|------------------|------------------|------------------|------------------|------------------|
|  | Thermal gradient      | Solvent gradient | Thermal gradient | Solvent gradient | Thermal gradient | Solvent gradient |
| Acetonitrile:Water                                 | 44.8                  | 39.7             | 15.8             | 15.2             | 32.7             | 32.5             |
| Acetonitrile:Water<br>with 0.1% TFA                | 45.4                  | 43.7             | 22.1             | 20.8             | 31.7             | 27.9             |
| Acetonitrile:<br>20 mM Ammonium<br>Hydroxide pH 10 | -                     | -                | 15.0             | 14.6             | 32.4             | 36.4             |


Comparable peak capacities were obtained for the solvent and thermal gradient runs.



# Efficiency Comparison of Solvent and Temperature Gradient for PRP<sup>®</sup>-1 Column




# Efficiency Comparison of Solvent and Temperature Gradient for Hypercarb<sup>®</sup> Column



◆ ◆ ◆ thermal gradient

◆ solvent gradient

# Efficiency Comparison of Solvent and Temperature Gradient for Blaze Column



◆ ◆ ◆ thermal gradient

◆ solvent gradient



# Comparison of Retention and Efficiency of PRP<sup>®</sup>-1 During High Temperature Evaluation

| Analyte                              | RT   | plates | area      |
|--------------------------------------|------|--------|-----------|
| uracil initial                       | 0.84 | 1,234  | 2,916,552 |
| uracil after initial runs before TFA | 0.86 | 1,300  | 3,061,639 |
| uracil after TFA before ammonia      | 0.80 | 1,357  | 3,118,760 |
| uracil after ammonia                 | 0.81 | 1,603  | 2,931,591 |
| phenol initial                       | 1.84 | 1,250  | 193,020   |
| phenol after initial runs before TFA | 1.87 | 1,263  | 190,543   |
| phenol after TFA before ammonia      | 1.82 | 1,230  | 203,936   |
| phenol after ammonia                 | 1.83 | 1,358  | 202,660   |



# Comparison of Retention and Efficiency of Hypercarb<sup>®</sup> During High Temperature Evaluation

| Analyte                              | RT   | plates | area      |
|--------------------------------------|------|--------|-----------|
| uracil initial                       | 1.43 | 3,827  | 3,152,477 |
| uracil after initial runs before NH3 | 1.45 | 5,227  | 2,698,373 |
| uracil after NH3 before TFA          | 1.45 | 4,658  | 2,991,397 |
| uracil after TFA                     | 1.47 | 5,298  | 2,909,248 |
| phenol initial                       | 2.41 | 5,414  | 182,970   |
| phenol after initial runs before NH3 | 2.56 | 5,466  | 213,160   |
| phenol after NH3 before TFA          | 2.36 | 4,811  | 209,921   |
| phenol after TFA                     | 2.63 | 4,944  | 207,484   |



# Comparison of Retention and Efficiency of Blaze C<sub>8</sub> During High Temperature Evaluation

| Analyte                   | RT   | plates | area      |
|---------------------------|------|--------|-----------|
| uracil initial            | 1.05 | 3,898  | 2,637,763 |
| uracil after initial runs | 1.06 | 4,826  | 3,112,145 |
| uracil after TFA          | 1.06 | 5,006  | 3,071,317 |
| phenol initial            | 1.81 | 7,162  | 166,090   |
| phenol after initial runs | 1.89 | 7,133  | 207,308   |
| phenol after TFA          | 1.89 | 7,591  | 201,861   |



# Conclusions

- Column bleed can be an issue with zirconia columns during temperature programming
- Several columns do give good performance under high temperature conditions without evidence of stationary phase degradation
- Similar retention characteristics are observed with solvent gradients and temperature programming
- Comparable or better peak quality is observed with temperature programming when compared to solvent gradients for the columns evaluated



# Acknowledgements

- **Thermo Hypersil-Keystone**
  - Steve Kozel and Rick Ludwig
- **Hamilton Company**
  - Dan Lee and Mike Benning





# Turn up the Heat!



**Selerity Technologies Inc.**

**2484 W. Custer Rd.**

**Salt Lake City, UT 84121**

**801-978-2295**

**[hotlc@selerity.com](mailto:hotlc@selerity.com)**

**[www.selerity.com](http://www.selerity.com)**

**Booth #208**

