

Natural Product Analysis Using Temperature Programmed HPLC: Sub-Ambient through 200 °C

Dale Felix, Stephanie J. Marin
Brian Jones, Jody Clark

Selerity Technologies, Inc.
Salt Lake City, UT
www.selerity.com



Temperature - the overlooked optimization parameter in HPLC

•Although nearly all of the physical parameters that play a role in liquid chromatographic separation are a function of temperature, temperature has not yet been adequately explored as a parameter to tune separation and shorten analysis times in LC.”*



*Nebojsa M. Djordjevic, Patrick W.J. Fowler, Fabrice Houdiere *J. Microcolumn Separations* 11(6) (1999) 403-413

The Use of Temperature in HPLC

- **Temperature Programming**
- **Speed**
- **Efficiency**
- **Decreased organic solvent usage**
 - ✓ Use less organic in solvent ratio
 - ✓ Perform isocratic separations and recycle solvent



Temperature Programming

- Temperature Programming in HTLC provides added benefits
 - Water becomes more like methanol as temperature increases so less organic solvent is required
 - Even faster analysis times*
 - Improves detectability for highly retained solutes*
 - Offers an alternative to gradient elution*



*Jian Chen, Dirk Steenackers, Andrei Medvedovici, Pat Sandra J.
High Res. Chromatogr. Vol 16 (1993) 605-608

Faster and More Efficient Separations

- **Speed**
 - Flatter van Deemter curves allow operation at flow rates many times optimal velocity
- **Higher efficiency - better resolution**
 - Increased diffusion rates provide lower plate heights at higher flow rates
 - Lower viscosity and back pressure permits higher flow rates with smaller particle size packings



Change in Retention with Temperature

$$k_1' / k_2' = \exp(\Delta H(T_2 - T_1) / (RT_1T_2))$$

Effect of T on retention factor depends on enthalpy of the solute. The larger the enthalpy, the greater the change in k'




Temperature and Elution Strength

- Increasing temperature 4 to 5 °C is comparable to increasing the methanol or acetonitrile concentration by 1% in a reversed phase system
- Viscosity is reduced 1 to 2 % per °C increase
 - Result – a significant reduction in back pressure

$$P = \eta \mu L / d_p^2$$




Solvent Polarity as a Function of Temperature



Data from Y. Yang et al. *J. Chromatogr. A* **810** (1998) 149.



Temperature Effects on Plate Height



Column Efficiency

- Efficiency increases with increasing temperature
 - Increased diffusivity with increasing temperature
 - Reduction of C term in vanDeemter equation




Polaratherm™ Design



- Temperatures from sub-zero to 200 °C
- Forced air circulation
- Isothermal and thermal gradient operation
- Integrated solvent preheating
- Effluent temperature control
- Flammable vapor sensor
- Compatible with any HPLC system

Separation of Barbiturates

Mobile Phase Preheating Improves Chromatography



Zirchrom PBD, 80°C


Fast temperature programming

- Obtain unique selectivities
- Replace solvent gradients with temperature gradients
- Simplifies the method development process



Natural Products

Ginseng Tablet Extract



Column: Selerity Blaze[®] C₈, 3 μm, 100 x 4.6 mm

Mobile Phase: 80:20 Acidic water:MeOH


Flow Rate: 0.70 mL/min

Detector: UV at 254

Temperature Program: 40 to 90 °C at 10°/min



Natural Products - Cocoa Extract



Column: Selerity Blaze[®] C₈, 3 μm, 100 x 4.6 mm

Mobile Phase: 80:20 Acidic water:MeOH


Flow Rate: 0.70 mL/min

Detector: UV at 254

Temperature Program: 40 to 90 °C at 10°/min



Natural Products - Green Tea Extract



Column: Selerity Blaze[®] C₈, 3 μm, 100 x 4.6 mm

Mobile Phase: 80:20 Acidic water:MeOH

Flow Rate: 0.70 mL/min

Detector: UV at 254


Temperature Program: 40 to 90 °C at 10°/min



Separation of Spotted Knapweed Extract

Column: Blaze C8,
150 x 4.6 mm **Mobile**
Phase: 20:80
ACN/Water
with 0.1%TFA
Flow Rate: 1.0 mL/min
Detection: UV 254 nm
Temperature Program:
hold at 40°C 2 min,
ramp to 100 at 30°/min,
hold five min


Sample Prep:
5 g weed extracted with
100 mL water (100 °C)




Separation of Russian Knapweed

Column: Blaze C8,
150 x 4.6 mm **Mobile**
Phase: 20:80
ACN/Water
with 0.1%TFA
Flow Rate: 1.0 mL/min
Detection: UV 254 nm
Temperature Program:
hold at 40°C 2 min,
ramp to 100 at 30°/min,
hold five min

Sample Prep:
5 g weed extracted with
100 mL water (100 °C)



Separation of Carotenoid Mixture at 15°C



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm
Mobile Phase: 1.25% THF in acetonitrile
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Temperature: 15°C, isothermal


70 minute run!

Elution Order:

Lutein
Zeaxanthin
Lycopene
β-carotene



Separation of Carotenoid Mixture at 25°C




Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm
Mobile Phase: 1.25% THF in acetonitrile
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Temperature: 25°C, isothermal

Elution Order:
Lutein
Zeaxanthin
Lycopene
 β -carotene




Separation of Carotenoid Mixture at 35°C



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm
Mobile Phase: 1.25% THF in acetonitrile
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Temperature: 35°C, isothermal
Elution Order:
Lutein
Zeaxanthin
Lycopene
 β -carotene



Separation of Carotenoid Mixture By Temperature Programming



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm

Mobile Phase: 1.25% THF in acetonitrile

Flow Rate: 1.0 mL/min

Detection: UV 450 nm

Temperature Program: hold at 15°C for five minutes, ramp to 35°C over two minutes, hold 15 minutes.



Conclusions

- Thermal gradients can replace solvent gradients for analysis of many natural products
- Increasing the temperature reduces the amount of organic modifier needed
- Temperature programming greatly simplifies method development for complex natural product separations




Acknowledgements


W. Jeffrey Hurst, Hershey Foods

John Estes and Neal Craft, Craft Technologies

Steven Burningham, Utah State Dept of Agriculture



Turn up the Heat!



Selerity Technologies Inc.

www.selerity.com



Booth #2204