

# **Temperature vs Pressure - Achieving High Resolution, Highly Efficient HPLC Separations Using Elevated Temperature In Place of High Pressure**

**Stephanie J. Marin, Dale Felix, Brian Jones and Jody Clark**

Selerity Technologies  
Salt Lake City, UT

**Pat Sandra, Frank David**

Research Institute for Chromatography  
Kortrijk, Belgium



# The Use of Temperature in HPLC

Temperature is considered to be the overlooked or forgotten optimization parameter in HPLC by many of the experts.

***“Although nearly all of the physical parameters that play a role in liquid chromatographic separation are a function of temperature, temperature has not yet been adequately explored as a parameter to tune separation and shorten analysis times in LC .”\****

\* Nebojsa M. Djordjevic, Patrick W.J. Fowler, Fabrice Houdiere *J. Microcolumn Separations* 11(6) (1999) 403-413



# Benefits of High Temperature Increased Diffusivity

- Increasing the temperature increases the enthalpy of solute transfer from mobile phase to stationary phase\*
  - Improves efficiency, particularly for large analytes
  - Allows operation at higher flow rates without penalty

\*F.D. Antia and Cs. Horvath, *J. Chromatogr.* **435** (1988) 1-15.

\*B. Yan, J. Zhao, J.S. Brown, J. Blackwell, P. W. Carr, *Anal. Chem.* **72** (2000) 1253-1262



# Benefits of High Temperature Decreased Viscosity

- As the temperature increases the viscosity of the eluent decreases thus lowering the system back pressure
  - Perform analysis at higher flow rates without over-pressurizing the pump
  - Use smaller particles to increase efficiency
  - Extend column length to add plates




# Benefits of Temperature Programming

- ❑ Allows increased plate count by adding column length
- ❑ Increases speed by increasing flow rate as viscosity decreases at elevated temperatures
- ❑ Changes selectivity as function of temperature and column phase
- ❑ Reduces organic solvent content – “Green Chemistry”




# Selerity Polaratherm™ Series 9000 Total Temperature Controller

- Forced air oven and chiller
- Isothermal and thermal gradient operation
  - Sub-zero to 200°C
  - Thermal gradients up to 30°C/min
- Mobile phase preheating and pre-cooling
- Peltier effluent temperature control
- Vapor sensor
- Compatible with any HPLC system
- Integrated software control
  - Waters, Agilent, and EZChrom




# High efficiency analyses

Utilization of multiple columns to improve efficiency



# Increase Resolution by Increasing Number of columns



**Column:** Agilent SB-C18 150 x 3.0 mm, 3.5  $\mu$ m

**Mobile phase:** ACN:Water 40:60

**Flow Rate:** 0.4 mL/min

**Temperature = 80°C**





# Efficiency for analytes using one, four and eight columns

Analyte	1 column	4 columns	8 columns
caffeine	5,000	41,000	46,000
pyridine	6,000	35,000	56,000
phenol	6,000	43,000	68,000
aniline	6,000	43,000	78,000
benzene	10,000	43,000	79,000
toluene	10,000	42,000	81,000

**Column:** Agilent SB-C18 150 x 3.0 mm, 3.5  $\mu$ m


**Mobile phase:** ACN:Water 40:60

**Flow Rate:** 0.4 mL/min

**Temperature = 80°C**



# Plate numbers vs columns




# Analysis of PAHs on ODS column

□ Column: **15 cm**

Zorbax SB-C18 150 x 3.0 mm, 3.5  $\mu\text{m}$   $d_p$

VWD1 A, Wavelength=254 nm (PAH-WIM\TEMP0013.D)

- 0.4 mL/min
- 80% to 100% methanol in 20 min
- 90°C




# Analysis of PAHs using four ODS Columns

□ Column: 60 cm

**4x** Zorbax SB-C18 150 x 3.0 mm, 3.5  $\mu\text{m}$   $d_p$


VWD1 A, Wavelength=254 nm (PAH-WIMKOPPEL04.D)




# Analysis of PAHs

## Aligned Chromatograms

VWD1 A, Wavelength=254 nm (PAH-WIM\KOPPEL04.D)





VWD1 A, Wavelength=254 nm (PAH-WIM\TEMP0013.D)




# Analysis of PAHs

## Resolution Comparison

4 Identical columns:  $N \times 4 \Rightarrow R_s \times 2$



# Aldehydes and Ketones of CARB 1004 using a solvent gradient



**Column:** Selerity Blaze<sub>200</sub> C<sub>18</sub>, 3 mm

150 x 4.6 mm columns

**Mobile Phase:** 10:90 acetonitrile:water, 10 to 60% ACN over 30 minutes.

**Flow Rate:** 1.5 mL/min

**Detection:** UV-VIS 365 nm

**Elution Order:**

Formaldehyde-DNPH

Acetaldehyde-DNPH

Acetone-DNPH

Acrolein-DNPH

Propionaldehyde-DNPH

Crotonaldehyde-DNPH

MEK-DNPH

Methacrolein-DNPH

Butyraldehyde-DNPH

Benzaldehyde-DNPH


Valeraldehyde-DNPH

Tolualdehyde-DNPH

Hexanaldehyde-DNPH



# Aldehydes and Ketones of CARB 1004 at 150°C



**Column:** Selerity Blaze<sub>200</sub> C<sub>18</sub>, 3 mm

2 150 x 4.6 mm columns

**Mobile Phase:** 25:75 acetonitrile:water

**Flow Rate:** 3.0 mL/min

**Detection:** UV-VIS 365 nm

**Temperature:** 150°C isothermal

**Elution Order:**

Formaldehyde-DNPH

Acetaldehyde-DNPH

Acetone-DNPH

Acrolein-DNPH

Propionaldehyde-DNPH

Crotonaldehyde-DNPH

Methacrolein-DNPH

MEK-DNPH

Butyraldehyde-DNPH

Benzaldehyde-DNPH

Valeraldehyde-DNPH


Tolualdehyde-DNPH

Hexanaldehyde-DNPH





# Aldehydes and Ketones of CARB 1004 at 150°C



**Column:** Selerity Blaze<sub>200</sub> C<sub>18</sub>, 3 mm

**2 150 x 4.6 mm columns**

**Mobile Phase:** 25:75 acetonitrile:water

**Flow Rate:** 3.0 mL/min

**Detection:** UV-VIS 365 nm

**Temperature:** 150°C isothermal

**Column:** Selerity Blaze<sub>200</sub> C<sub>18</sub>, 3 mm

**2 150 x 4.6 mm columns**

**Mobile Phase:** 10:90 acetonitrile:water, 10 to 60% ACN over 30 minutes.

**Flow Rate:** 1.5 mL/min

**Detection:** UV-VIS 365 nm



# Conclusions

- ❑ High resolution separations can be achieved by increasing temperature.
- ❑ Lower mobile phase viscosity and system back pressure allow the use of longer columns to increase efficiency and resolution at elevated temperature.

# Acknowledgements

Dr. Jeffrey Loo, GM



# Turn up the Heat!



Stop by booth **6452**  
or visit our website  
**[www.selerity.com](http://www.selerity.com)**

