

The Effect of Temperature on Selectivity in HPLC

Brian Jones, Jody Clark, Dale Felix and Stephanie J. Marin

Selerity Technologies, Inc.

2484 W. Custer Road

Salt Lake City, UT 84104

www.selerity.com



The Use of Temperature in HPLC

Temperature is considered to be the overlooked or forgotten optimization parameter in HPLC by many of the experts.

“Although nearly all of the physical parameters that play a role in liquid chromatographic separation are a function of temperature, temperature has not yet been adequately explored as a parameter to tune separation and shorten analysis times in LC .”*



* Nebojsa M. Djordjevic, Patrick W.J. Fowler, Fabrice Houdiere *J. Microcolumn Separations* 11(6) (1999) 403-413

Typical Retentive Behavior

- Retention uniformly decreases with temperature
- Linear van't Hoff plots – $\log k$ vs. $1/T$
- A 1% methanol increase is equivalent to $\sim 4^{\circ}\text{C}$ temperature increase



Change in Retention with Temperature

$$k_1' / k_2' = \exp(\Delta H(T_2 - T_1) / (RT_1T_2))$$

- Effect of temperature on retention factor depends on enthalpy of the solute. The larger the enthalpy, the greater the change in k'



Analyte Families Giving Reversed Retention with Increasing Temperature

- Polyethers
 - Reverse retention due to apparent decreasing polarity due to hydrogen bond strength reduction
- Some dipeptides
 - Decreasing polarity with conformational change




Factors Affecting Retention That Are Influenced by Temperature

- Hydrogen bonding
- Solvation sphere around analytes
- Hydration extent of column surface
- Functional group interaction
- Ordering, shape and hydration transitions
 - Conformation changes
- Dielectric constant of mobile phase




Temperature Affects Hydrogen Bonding

- Increasing temperature-
 - Increases intermolecular distance
 - Weakens hydrogen bonds




Column Surface Hydration




Functional Group Conformation

Octadecyl ordering at low temperatures




Functional Group Conformation

Octadecyl disorder at higher temperatures



Functional Group Phase Transitions

- *N*-Isopropylacrylamide coated silica




Analyte Conformational Changes

- Heteroduplex and homoduplex oligonucleotides separable by HPLC with temperature gradients
 - Axial or spatial gradients
 - Temperature programming of the entire column

From R.E. Gerber and R.G. Hatch in US Patent 6,486,309



Solvent Polarity as a Function of Temperature



High Temperature Liquid Chromatography Advantages

- The advantages of using temperature to optimize separations in HPLC are well documented in the literature.
- The major advantages are
 - Increased speed
 - Higher efficiencies and resolution
 - Ability to tune selectivity with temperature
 - Decreased organic solvent usage
 - Use less organic in solvent ratio
 - Perform isocratic separations and recycle solvent



Increased Diffusivity

- Increasing the temperature increases the enthalpy of solute transfer from mobile phase to stationary phase*
 - Improves efficiency, particularly for large analytes
 - Allows operation at higher flow rates without penalty



*F.D. Antia and Cs. Horvath, *J. Chromatogr.* **435** (1988) 1-15.

*B. Yan, J. Zhao, J.S. Brown, J. Blackwell, P. W. Carr, *Anal. Chem.* **72** (2000) 1253-1262

Decreased Viscosity

- As the temperature increases the viscosity of the eluent decreases thus lowering the system back pressure
 - Perform analysis at higher flow rates without over-pressurizing the pump
 - Use smaller size packing materials in columns increasing efficiency




Temperature and Elution Strength

- Increasing temperature 4 to 5 °C is comparable to increasing the methanol or acetonitrile concentration by 1% in a reversed phase system
- Viscosity is reduced 1 to 2 % per °C increase
 - Result – a significant reduction in back pressure

$$P = \eta \mu L / d_p^2$$



Elution Strength Illustration: Separation of Steroids Using Water as the Mobile Phase



Column: ZirChrom PBD, 3 μ m
100 X 4.6 mm

Detection: UV 254 nm

Flow Rate: 6.0 mL/min

Mobile Phase: Water

Temperature: 200°C

Flow Rate: 3.0 mL/min

Mobile Phase: 25:75 acetonitrile:water

Temperature: 50°C

Elution Order:

Uracil


Androstadienedione

Androstenedione

Epitestosterone




Prototype C₁₈ Silica Column Separation of Steroids



Efficiency with High Speed

- Temperature effect on plate height




Historical Obstacles in the Pathway of Extended Range HPLC

- Most commonly listed reasons why temperature has not been utilized as an optimization tool in reversed phase
 - Poor temperature stability of silica based columns
 - Lack of adequate column heating system
 - Mobile phase composition is easily adjusted and provides the flexibility to handle a wide range of samples



Mobile Phase Pre-heating




Stable Column Phases

- Silica (Selerity Blaze C₈)
- Zirconia* (ZirChrom PBD, DiamondBond, CARB)
- Carbon (Thermo Hypercarb)
- Polymeric (Jordi, Hamilton with SS fittings)

* Not recommended for temperature programmed conditions because of severe column bleed



Blank Runs with Thermal Gradients



50:50 Acetonitrile:Water at 254 nm



Polaratherm™ Design

- Temperatures from sub-zero to 200 °C
- Forced air circulation
- Isothermal and thermal gradient operation
- Integrated solvent preheating
- Effluent temperature control
- Flammable vapor sensor
- Compatible with any HPLC system



Diesel Range Aromatics - Columns

- **Partisil Amino-Cyano (PAC) column**
 - **5 μ m, 250 x 4.6 mm**
- **All stainless steel hardware**
- **Hexane mobile phase**




Retention Data for PAC

Column Isothermal Conditions

Analyte	Retention Time			
	0°C	30°C	60°C	100°C
toluene	4.56	4.23	3.89	3.54
tetralin	4.59	4.29	4.00	3.69
thiophene	5.02	4.46	4.08	3.69
naphthalene	6.81	5.69	4.99	4.32
acenaphthene	7.49	5.89	4.92	4.47
benzothiophene	7.81	5.82	5.22	4.49
dibenzothiophene	11.25	8.38	6.73	5.47
anthracene	12.14	9.04	7.30	5.86



Retention Data for PAC Column Isothermal Conditions



β -Carotene and Isomers at 24 °C




ACN/THF 90:10
Isocratic
1 mL/min
Column: Genesis C₁₈




β -Carotene and Isomers at 10 °C



β -Carotene and Isomers at 0 °C




Separation of Carotenoid Mixture at 15°C



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm
Mobile Phase: 1.25% THF in acetonitrile
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Temperature: 15°C, isothermal
70 minute run time



Separation of Carotenoid Mixture at 25°C



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm

Mobile Phase: 1.25% THF in acetonitrile


Flow Rate: 1.0 mL/min

Detection: UV 450 nm

Temperature: 25°C, isothermal




Separation of Carotenoid Mixture at 35°C



Column: Jones
Chromatography Genesis
C₁₈, 150 x 4.6 mm
Mobile Phase: 1.25%
THF in acetonitrile
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Temperature: 35°C,
isothermal



Separation of Carotenoid Mixture Using a Temperature Program



Column: Jones Chromatography Genesis C₁₈, 150 x 4.6 mm

Mobile Phase: 1.25% THF in acetonitrile


Flow Rate: 1.0 mL/min

Detection: UV 450 nm

Temperature Program: hold at 15°C for five minutes, ramp to 35°C over two minutes, hold 15 minutes.



Separation of Analgesics on a Selerity Blaze™ C₈ Using a Thermal Gradient



Column: Selerity Blaze C₈, 3 μm
100 x 4.6 mm

Mobile Phase: 40:60 acetonitrile:water
with 0.1%TFA

Flow Rate: 1.5 mL/min

Detection: UV 220 nm

Temperature Program: hold at 50°C for
one minute, ramp to 100°C at 30°C/min,
hold six min.

Elution Order:

Acetaminophen

Caffeine

Salicylamide

Aspirin


Salicylic acid

Ibuprofen

Naproxen



Analgesics Using a Hypercarb[®] Column and a Thermal Gradient



Column: Thermo Electron Hypercarb[®], 7 μ m, 100 x 4.6 mm
Mobile Phase: 35:65 acetonitrile:water with 0.1% TFA
Flow Rate: 4.0 mL/min
Detection: UV 220 nm
Temperature Program: thermal gradient from 125° to 200°C at 30°/min, hold five min.

Elution Order:

Caffeine
Aspirin
Salicylic Acid
Ibuprofen
Phenacetin
Acetaminophen
Naproxen



Conclusions

- Temperature is powerful but under-used tool for selectivity tuning.
- Increasing the temperature reduces the amount of organic modifier needed and generally speeds elution.
- Temperature adjustment can be a fast way to optimize a separation.
- Shorter analysis times with better efficiency usually result from higher temperatures.



Acknowledgements

The authors thank Thermo Electron and Hamilton Company for providing columns for this work.

Appreciation is also given to Craft Technologies for the carotenoid data.



Turn up the Heat!



Selerity Technologies Booth 2204