

Use and Practicality of Water-Only HPLC Separations

Jody Clark, Stephanie J. Marin, Brian A. Jones and
W. Dale Felix

Selerity Technologies, Inc.
2484 W. Custer Road
Salt Lake City, UT 84104



www.selerity.com

Introduction

- Properties of water
- Advantages of using temperature
- Rules of thumb for determining the feasibility of using pure water as the mobile phase for separations
- Practical applications




Why Use Water as a Mobile Phase?

- Inexpensive
- Readily available
- Non-polluting
- Transparent to most detectors including UV, FID, and NMR (D_2O)
- Elution strength can be controlled by varying temperature



Unique Properties of Water are Primarily Due to Hydrogen Bonding

- Increasing temperature
 - Increases intermolecular distance
 - Weakens hydrogen bonds



- Results in
 - Decreasing density and viscosity



Structure of High Temperature Water (HTW)

- HTW is structurally different from ambient liquid water
 - An infinite percolating network of H-bonds exist in ambient water
 - Small clusters of H-bonded water molecules exist in HTW
 - As the temperature increases and the density decreases the average cluster size decreases
 - The breaking of the H-bond network reduces the barrier for translational and rotational motion




Changes in the Dielectric Constant

- Changes in the extent of hydrogen bonding are accompanied by corresponding changes in the dielectric constant
 - With increasing temperature and decreasing density, the dielectric constant of water decreases
 - HTW behaves more like polar organic solvents
 - Small organic compounds are highly soluble in HTW




Solvent Polarity as a Function of Temperature



Data from Y. Yang et al. *J. Chromatogr. A* **810** (1998) 149.



Viscosity of Water vs. Temperature




High Temperature Liquid Chromatography (HTLC) Advantages

- Decreasing analysis time
 - Due to ease of increasing linear velocity
 - Due to a generally decreasing retention for most compounds
- Increasing efficiencies and resolution
- Selectivity tuning
- Decreasing organic solvent usage



Chromatographic Example of Using High Temperature



Column: ZirChrom PBD, 3 μ m
100 X 4.6 mm

Detection: UV 254 nm

Flow Rate: 6.0 mL/min

Mobile Phase: Water

Temperature: 200°C

Flow Rate: 3.0 mL/min

Mobile Phase: 25:75 acetonitrile:water

Temperature: 50°C

Elution Order:

Uracil


Androstadienedione

Androstenedione

Epitestosterone



Instrument Diagram



Columns Available for Water Separations

- Polymer (temp. limit: 150-175°C)
- Graphitic Carbon (temp. limit: above 200°C)
- Bridged Ethyl Hybrid (temp. limit: 200°C)
- Selected Silica-highly cross linked (temp. limit: 100-150°C)




Rules of Thumb


- Analytes that are good candidates for using water as the mobile phase
 - Soluble in water
 - Slightly soluble in water
 - Insoluble in water but soluble in alcohols, glycerin, and acidified water



Glycols on 1.0mm ID Hypercarb™



Analytes Soluble in Water (found in diet cola)



Column: Xbridge C8, 3.5 μm
150 X 2.1 mm

Detection: UV 214 nm

Flow Rate: 1.5 mL/min

Mobile Phase: Water

Temperature: 110°C

Elution Order:

Aspartame


Caffeine

Sodium Benzoate

Acesulfame K



Addition of ~0.005% Formic Acid and Temperature Programming



Column: Xbridge C8, 3.5 μ m
150 X 2.1 mm

Detection: UV 214 nm

Flow Rate: 1.5 mL/min

Mobile Phase: 0.005% formic acid in water

Temperature: 95°C

hold 2 min.

ramp to 110°C @ 10°C/min.

Elution Order:

Aspartame





Caffeine

Sodium Benzoate

Acesulfame K




Compounds Slightly Soluble

Name	Structure	Solubility
Pyrogallol		~590 g/L
Catechol		~450 g/L
Phenol		~80 g/L
Guaiacol		~16 g/L



Hydroxyphenols Eluting in Order of Analyte Affinity for Water



Column: Polymer Labs PLRP-S, 3.0 μm
150 X 2.1 mm

Detection: UV 254 nm

Flow Rate: 1.5 mL/min

Mobile Phase: Water

Temperature: 175°C

Elution Order:

Pyrogallol


Catechol

Phenol

Guaiacol



Alcohol Based Mouthwash



Column: Xbridge C18, 3.5 μm
150 X 2.1 mm

Detection: UV 254 nm
FID @ 400°C

Flow Rate: 1.0 mL/min

Mobile Phase: Water

Temperature: 200°C

Split Ratio: 9:1 with 0.1 mL/min into FID

Elution Order:

Methyl Salicylate


Thymol

Eucalyptol

Menthol



Flexibility with Temperature Programming (Alcohol Based Mouthwash)



Column: Xbridge C18, 3.5 μ m
150 X 2.1 mm

Detection: UV 254 nm
FID @ 400°C

Flow Rate: 1.5 mL/min

Mobile Phase: Water

Temperature: 125°C

hold 1 min.

ramp to 200°C @ 25°C/min.

Split Ratio: 14:1 with 0.1 mL/min into FID

Elution Order:

Methyl Salicylate


Thymol

Eucalyptol

Menthol



Glycerin Based Mouthwash



Column: Xbridge C18, 3.5 μ m
150 X 2.1 mm

Detection: UV 254 nm
FID @ 400°C

Flow Rate: 1.5 mL/min

Mobile Phase: Water

Temperature: 135°C

hold 1 min.

ramp to 200°C @ 25°C/min.

Split Ratio: 14:1 with 0.1 mL/min into FID



Most Practical Use of Using Water as the Mobile Phase

- Gain use of the flame ionization detector (FID)
 - Mass sensitive detector
 - Response is proportional to the number of carbon atoms being burned
 - Wide linear range, about 10^8
 - Responds with high sensitivity to organic compounds



Conclusion

- The use of water as a mobile phase is possible when used with temperature for compounds that are soluble in water, alcohol and glycerin
- Great choice when analyzing compounds that do not contain chromophores



Acknowledgements

- Nelu Grinberg, Ph.D., Boehringer Ingelheim Pharmaceuticals, Inc.

